CROWN & BRIDGE



- 900 mpa strength
- High translucency
- Chamfer prep with RILA and 2mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 2 units
- Polished after milling to remove micro ridges
- Lifetime warranty

≅e.max[°]

PRESSED ANTERIOR

- 450 mpa strength
- 97% clinical rating
- 96.8% survival rate

MILLED POSTERIOR

- 500 mpa strength
- 97% clinical rating
- 96.8% survival rate



- >1200 mpa strength
- Medium translucency
- Light Chamfer prep with 1.5mm occlusal clearance
- Intaglio groove for better retention
- Pontic span up to 4 units
- Polished after milling to remove micro ridges
- Lifetime warranty



IVOCLAR CLASSIC PORCELAIN TO ARGEN METALS (PFM)

<u>2202LifeUSA</u> + <u>**≊e.max**</u>

PORCELAIN FUSED TO ZIRCONIA (PFZ)

FULL - CAST RESTORATIONS



5968 NC-16 Business Hwy | Maiden | NC | 28650 keydentalnc.com | 704.489.9978